Stahlzement

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Steel cement

04090	Powder light yellow, jar	40 g
04091	Powder yellowish white, jar	40 g
04093	Powder light grey, jar	40 g
04132	Liquid normal, bottle	20 ml
04133	Liquid fast, bottle	20 ml

Stone cement

04080	Powder light yellow, jar	40 g
04081	Powder yellowish white, jar	40 g
04083	Powder light grey, jar	40 g
04130	Liquid normal, bottle	20 ml
04131	Liquid fast, bottle	20 ml

DETAX

GmbH & Co. KG

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

Telefon 07243/510-0 · Fax 07243/510-100

www.detax.de · post@detax.de

Made in
Germany

03/2012

Drala® silico-phosphate cements

GB

Instructions for use

Stone Cement · Steel Cement

Dental water-based cements

Silico-phosphate cements

Filling cements

Mixing ratio

1,5 g of powder + **0,4 ml** of liquid

The recommended by-weight ratio is 2,4 parts of powder to 1 part of liquid.

Mixing

In order to achieve correct mixing of powder and liquid, a glass slab with a temperature of 18° - 21°C / 64° - 70 °F and a plastic or agate spatula should be used. Be careful to use only clean dry glass slab and spatula. The cement should not come into contact with any moisture! When mixing the cement, first a little of powder is added to the liquid and then gradually more until the correct consistency is reached. The powder is mixed in by using alternately circular and patting movements. Fresh amounts of liquid should never be added during mixing!

The filling consistency is reached when the paste seems to be rather dry and, on lifting the spatula, tears off.

Please note

- Increased plasticity or a longer setting period should not be obtained by using a thinner consistency but by cooling the mixing slab. A too thin mixture may impair the quality of the cement.
- Overdosing of powder or increased temperatures accelerate the setting time.

Application

The filling should be realized immediately after mixing, by means of a sili-coresistant instrument. Be sure to avoid penetration of moisture!

Caution

The liquid contains phosphoric acid that may cause acid burns. Avoid contact with eyes, mucosa and skin. In case of accidental contact immediately flush affected body areas with generous amounts of water and soap. In case of contact with eyes immediately flush under running water for at least 15 minutes and seek immediately medical care.

Storage

Close recipients carefully after each use. Keep tightly closed at temperatures of 15 °C - 25 °C/59 °F - 77 °F.

Keep away from direct sunlight.

Only for intended use by qualified personnel.

Indications

- filling cement for molars and premolars, if a sufficient pulp protection is existing.

Contra indications

- Drala® Silico-phosphate cement is not allowed to be applied directly in close proximity to the unprotected pulp.

Attention

When working in close proximity to the pulp a calcium hydroxyde lining should be applied to the cavity before applying the silico-phosphate cement. Also the cement should not be applied in case of humid working surface and for fillings under intense masticatory pressure.

Mixing time

- 1 minute

Working time

- 2 minutes*

Setting time

- normal setting:
4 - 4½ minutes**
- fast setting:
3 - 3½ minutes**

* from beginning of mixing, at 21 °C/70 °F, 50 ± 5 % rel. humidity

** from end of mixing, at 37 °C/99 °F, 50 ± 5 % rel. humidity

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